

Live life in high Res

Nasal and nasal pillows
range of masks



ResMed

Drift to sleep wearing our most comfortable gear

Let's find a mask that fits your needs

At ResMed, we believe a mask should seal securely and comfortably – no matter how you sleep or breathe. That's why each of our masks is designed with different features to meet a variety of needs. After all, every person and every face is unique. With the number of options we have available, we're confident you'll discover one that feels tailor-made.



NASAL
MASK

What is a nasal mask?

Most nasal masks consist of a cushion that covers the nose between the nasal bridge and upper lip. Some nasal masks have cushions that cradle under the nose. Both types are great options if you:

- Breathe through your nose
- Prefer less coverage on your face



MONEY BACK GUARANTEE

We're so sure that you will love our AirFit™ nasal and nasal pillows masks that we are offering a 30 day money back guarantee if you are not completely satisfied.*

*Terms and conditions apply. See in store or online for details.

Enjoy a relaxing night's sleep wearing the latest from our nasal and nasal pillows mask range.



NASAL
PILLOWS MASK

What is a nasal pillows mask?

A nasal pillows mask is a lightweight and compact mask with tubular silicone inserts that fit gently inside your nostrils to deliver therapy air. Nasal pillows masks are a great option for:

- Reading or watching TV in bed
- Reducing facial markings



Rest assured, we'll help you find the right mask

Our nasal and nasal pillows mask range will continue to grow and get better, and with our new mask categories, finding the right mask has never been easier. These categories make it easy for you to understand the key benefits of each mask, so you can take one home with confidence. Each category offers something a little different. Have a look to see which one suits you best.

Versatile Fit Masks

Designs that fit nearly every face.


Minimalist Masks

Small. Light. Simple.

Freedom Masks

Sleep and move comfortably.

Speak to your ResMed Authorised Dealer to find the right mask that will help you adjust to therapy, and to help achieve optimum results night after night.



AirTouch™ N30i

Freedom Fit

AirTouch N30i is ResMed's first tube-up nasal mask featuring a ComfiSoft™ fabric cushion and a frame that is fully wrapped in soft, breathable fabric.

AirTouch N30i is available in small and standard frame sizes. Both come with medium and small-wide cushions. Large cushion is available separately.

- 1 ComfiSoft™ nasal cradle cushion with QuietAir vent fuses fabric and silicone to create a soft and breathable cushion designed to enhance comfort.
- 2 SpringFit™ frame is fully wrapped in soft, breathable, moisture-wicking fabric to keep your skin feeling cool whilst also auto-adapting to different facial profiles.
- 3 Quick-release swivel elbow and no tubing in front means you can sleep in the position you want.





AirFit™ N30i

Nasal cradle mask

DID YOU KNOW?

This mask's frame and cushions are interchangeable with AirTouch N30i.


Freedom

ResMed AirFit N30i features a tube-up design that makes it easier to sleep in many positions and closer to your bed partner. It's a great option if you sleep on your side or stomach.

AirFit N30i nasal cradle cushions are interchangeable with AirFit P30i nasal pillows cushions. AirFit N30i is available in small and standard frame sizes. Both come with small, medium and small-wide cushions.

- 1 SpringFit frame adapts to different facial profiles to provide a personalised fit.
- 2 Soft nasal cradle cushion with QuietAir diffused venting to enhance comfort and provide a quieter night's sleep.
- 3 No front-end tube means the mask stays out of sight and out of the way.





AirFit™ P30i

Nasal pillows mask

DID YOU KNOW?

This mask's frame and cushions are interchangeable with AirTouch N30i.

Freedom

With a tube-up design and a nasal pillows cushion, AirFit P30i makes it easier for you to sleep in many different positions and with less mask on your face.

AirFit P30i nasal pillows cushions are interchangeable with AirFit N30i nasal cradle cushions. AirFit P30i is available in small and standard frame sizes. The small starter pack comes with small and medium cushions; the standard starter pack comes small, medium and large size cushions.

- 1 SpringFit frame adapts to different facial profiles to provide a personalised fit.
- 2 QuietAir vent is designed to minimise airflow and night-time disruptions.
- 3 No front-end tube means the mask stays out of sight and out of the way.



AirFit™ N30

Nasal cradle mask



Minimalist

ResMed's lightest mask ever¹, AirFit N30, is designed to make staying on therapy easy. Small and easy to use, it features a curved nasal cradle cushion and adjustable QuickFit™ headgear.

AirFit N30 includes small, medium and small-wide cushions.


- 1 Adjustable QuickFit headgear is designed to make fitting and adjustment easy.
- 2 Nasal cradle cushion curves under the nose for a comfortable and secure fit.
- 3 QuietAir™ vent technology diffuses exhaled airflow gently and quietly.



¹As of October 2019. The total weight of AirFit N30 is 45g, inclusive of the headgear but non-inclusive of the packaging.

AirFit™ P10

Nasal pillows mask



Minimalist

AirFit P10's sleek styling provides a real sense of facial freedom, and with just three parts, it's designed to be simple to use and maintain. AirFit P10 is our best-selling nasal pillows mask for several reasons.

AirFit P10 includes small, medium and large nasal pillows cushions, and it comes with a standard size headgear. It is also available in a for Her variant with smaller size option and uniquely styled headgear. Compatible with AirMini™.

- 1 Adjustable QuickFit headgear with the split-strap design is easy to put on, take off and adjust.
- 2 Nasal pillows with dual-wall technology seals comfortably on contact.
- 3 QuietAir™ woven-mesh venting diffuses air gently to reduce night-time disruptions.





AirFit™ N20

Nasal mask

Versatile Fit


AirFit N20 is designed to fit a broad range of face shapes and sizes. Thanks to its under-eye frame for greater visual freedom and all-round comfort, AirFit N20 has become our most popular nasal mask since its release.

AirFit N20 is available in medium and large sizes. It is also available in a for Her variant with smaller size option and uniquely styled headgear. Compatible with AirMini.

- 1 InfinitySeal™ silicone cushion is designed to prevent blow-outs and seal well under variable therapy pressures.
- 2 Soft frame padding has been incorporated into the face-side of the frame to help provide comfort at every touchpoint.
- 3 Magnetic clips guide the headgear to the frame quickly and easily.*



*The mask contains magnets that may interfere with certain implants or medical devices. Please refer to the User Guide for complete labelling information, including magnet contraindications and warnings.



AirTouch™ N20

Nasal mask

Versatile Fit

AirTouch N20 features an UltraSoft™ memory foam cushion, which is designed to provide a light, breathable seal. It shares the same frame and headgear as AirFit N20 so users can swap between cushions if desired.

AirTouch N20 is available in medium and large sizes. It is also available in a for Her variant with smaller size option and uniquely styled headgear.

- 1 UltraSoft memory foam is a softer alternative than silicone and designed to reduce the risk of red marks and skin irritation.
- 2 A disposable cushion (to be replaced monthly) helps make maintaining hygiene and cleanliness easier than ever before.
- 3 Magnetic clips guide the headgear to the frame quickly and easily.*



14
DAY

REPLACEMENT
GUARANTEE

If you're not completely satisfied with the UltraSoft memory foam cushion, you may exchange the cushions that are a part of the AirTouch N20 Starter Kit for an AirFit N20 InfinitySeal™ silicone cushion within 14 days of purchase. T&Cs apply. See in store or online for details.



BECOME A VIP TODAY...

Join sleepvantage and be treated to a range of exclusive member benefits.[^]

- ✓ **Members-only warranty****
Additional sleepvantage membership warranty on ResMed products.
- ✓ **Feel supported**
Access to our friendly product support line.
- ✓ **Stay informed**
Access to a range of educational articles and videos.
- ✓ **Get special offers**
Be part of exciting promotions and special offers.

For more information, contact sleepvantage today



sleepvantage.com.au



1300 305 705 (AUS) or 0800 737 633 (NZ)



info@sleepvantage.com.au

[^]sleepvantage membership is available to residents of Aus and NZ.

**Terms and conditions apply. See website for more information.



ResMed Pty Ltd
1 Elizabeth Macarthur Drive
Bella Vista NSW 2153 Australia

AirFit, AirMini, AirTouch, ComfiSoft, InfinitySeal, QuickFit, QuietAir, sleepvantage, SpringFit and UltraSoft are trademarks and/or registered trademarks of the ResMed family of companies. For patent information, see ResMed.com/ip. Specifications may change without notice. © 2022 ResMed Pty Ltd. 10111537r6 2024-10